

Dressing Application

1. Always use aseptic techniques when dressing a wound
2. Cut, fold or layer, the dressing to fit the wound
3. Loosely fill deep wounds, to allow the Kytocel to absorb wound fluid and to expand filling the wound. Do not pack tightly.
4. When using Kytocel ribbon dressings in deep cavity wounds, insert the ribbon in one piece and leave at least 2.5cm outside the wound for ease of removal
5. Cover and secure with an appropriate secondary dressing
6. For dry wounds (to aid a moist wound environment) place a Kytocel on the wound and wet with sterile water or saline over the wound area only. Cover the dressing with an occlusive or semi-occlusive dressing to maintain moisture at the wound site and to avoid subsequent dressing adherence to the wound
7. For minor bleeding wounds, apply the Kytocel to the wound with a light pressure, maintain the pressure for few minutes until the bleeding is stopped, secure the dressing with a suitable secondary dressing

Dressing Change and Removal

1. Kytocel shall be changed when the dressing has reached its absorption capacity, or when the conditions of the wound or the surrounding skin or the level of exudate requires.
2. If simple removal of the dressing is difficult due to adherence, it should be fully saturated with sterile saline or sterile water to fully gel the dressing, remove slowly
3. Irrigate the wound site with a suitable wound cleanser (e.g. sterile saline) prior to application of a new dressing

Storage

Store in a cool, dry place.

Do not re-sterilise

The use by date of this product is printed on the packaging.

SYMBOLS USED ON LABELLING



Single use



Do not Re-sterilize



Use by: year and month.



Sterile. Method of sterilisation: EO



Product Reference



Batch number.



Attention, see instructions for use.



Does not contain Latex



Storage Conditions



Manufacturer



EC Representative



CE mark



Keep Away from Sunlight



Keep Dry



Do not use if package is open or damaged



Medical Device



See instruction for use